

https://selectapi.datascope.refinitiv.com/RestApi/v1/Extractions/ExtractWithNotes

Save

POST https://selectapi.datascope.refinitiv.com/RestApi/v1/Extractions/ExtractWithNotes

Send

Params Authorization Headers (10) Body Pre-request Script Tests Settings

Cookies

Type Inherit auth fr...

The authorization header will be automatically generated when you send the request. Learn more about [authorization](#)

This request is not inheriting any authorization helper at the moment. Save it in a collection to use the parent's authorization helper.

Body Cookies Headers (16) Test Results

Status: 401 Unauthorized Time: 84 ms Size: 620 B Save Response

Pretty Raw Preview Visualize JSON

```
1 {
2   "error": {
3     "message": "The Authorization header is malformed. Supported Authorization types are Token,Basic."
4   }
5 }
```

https://selectapi.datascope.refinitiv.com/RestApi/v1/Extractions/ExtractWithNotes

POST https://selectapi.datascope.refinitiv.com/RestApi/v1/Extractions/ExtractWithNotes

Params Authorization Headers (10) Body Pre-request Script Tests Settings

<input checked="" type="checkbox"/>	Content-Length	<calculated when request is sent>	
<input checked="" type="checkbox"/>	Host	<calculated when request is sent>	
<input checked="" type="checkbox"/>	User-Agent	PostmanRuntime/7.29.2	
<input checked="" type="checkbox"/>	Accept	*/*	
<input checked="" type="checkbox"/>	Accept-Encoding	gzip, deflate, br	
<input checked="" type="checkbox"/>	Connection	keep-alive	
<input checked="" type="checkbox"/>	Prefer	respond-async	
<input checked="" type="checkbox"/>	Authorization	_B1izVzLKZHKgNcAEAqzvd1XYRrh7LgRebo7cymn3QhXP...	
	Key	Value	Description

Body Cookies Headers (16) Test Results

Status: 401 Unauthorized Time: 84 ms Size:

Pretty Raw Preview Visualize JSON

```
1 {
2   "error": {
3     "message": "The Authorization header is malformed. Supported Authorization types are Token,Basic."
4   }
5 }
```

https://selectapi.datascope.refinitiv.com/RestApi/v1/Extractions/ExtractWithNotes

Save

POST https://selectapi.datascope.refinitiv.com/RestApi/v1/Extractions/ExtractWithNotes

Params Authorization Headers (10) Body Pre-request Script Tests Settings

none form-data x-www-form-urlencoded raw binary GraphQL JSON

```
1  {}
2  ..... "ExtractionRequest": {}
3  ..... "@odata.type": "#DataScope.Select.Api.Extractions.ExtractionRequests.ElektronTimeseriesExtractionRequest",
4  ..... "ContentFieldNames": [
5  .....     "Ask",
6  .....     "Bid",
7  .....     "High",
8  .....     "Last",
9  .....     "Low",
10 .....     "Trade Date"
11 ..... ],
12 ..... "IdentifierList": {}
13 ..... "@odata.type": "#DataScope.Select.Api.Extractions.ExtractionRequests.InstrumentIdentifierList",
```

Body Cookies Headers (16) Test Results

Status: 401 Unauthorized Time: 80 ms Size: 614 E

Pretty Raw Preview Visualize JSON

```
1  {}
2  "error": {}
3  | "message": "The Authorization header is malformed. Supported Authorization types are Token,Basic."
4  | }
5  | }
```